vivocare
health enjoyment

Model: **Smilky**

EN Electric Breast Pump

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PT Bomba tira-leite elétrica

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DE Elektrische Muttermilchpumpe

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EN **Instruction Manual**

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Product Description	02
Safety Notice	03
Unit Illustration	05
Installation Instructions	06
Product Assembly	06
Battery Installation	10
The Power Adapter	11
Using the breast bump	12
Pumping	13
Size of Breastshield	15
Cleaning	16
Trouble Shooting	17
Waste Treatment	18
Maintenance	19
Specifications	21
Warranty	23
Electromagnetic Compatibility Information	24

Product Description

Thank you for purchase and use Smilky electric breast pump. The electric breast pump is electric powered suction device used to express and collect breast milk from lactating mother.

This product follows a baby's natural nursing rhythm. The unique design can help you to express breast milk quickly and comfortably. Please read this manual carefully before using stimulation and expression mode. Press the power button, the breast pump automatically starts up in stimulation mode. This is a single use product that should not be shared between moms.

Two mode technology of Stimulation and Expression :

This product has stimulation and expression mode. Press the power button, the breast pump automatically starts up in stimulation mode.

If no action is taken during stimulation mode, after two minutes, the pump will automatically change to the expression mode.



Stimulation mode: Fast sucking/pumping rhythm to stimulate the milk ejection reflex and to start the milk flowing.

Expression mode: Slower sucking/pumping rhythm for gentle and efficient milk removal as quickly as possible.

Intended use :

1. This breast pump is for household use only. As the pump is compact and discrete to use, you can take it with you anywhere, allowing you to express milk at your own convenience and maintain your milk supply.
2. If your breasts are engorged (painful or swollen), you can express a small amount of milk before or between feeds to ease the pain and to help your baby latch on more easily.
3. If you are separated from your baby and wish to continue to breastfeed when you are reunited, you should express your milk regularly to stimulate your milk supply.

Explanation of symbols:

SYMBOL EXPLANATION	
	Type BF applied part
	Caution
	Consult Accompanying Documents
	Keep Away From Sunlight
	Storage and Transportation Temperature Limit: -20°C~55°C(-4°F~131°F)
	Storage and Transportation Humidity Limit: 15~93%RH
	Serial Number
	European Authorized Representative
	Manufacturer
	Direct Current

Caution

1. Use the product only for its intended use as described in this manual.
2. Please check that all the parts of breast pump is in before using.
3. Do not use attachments not recommended by the manufacturer.
4. Do not attempt to remove the breastshield from your breast while pumping.
Turn the breastpump off and break the seal between your breast and breastshield with your finger, then remove breastshield from your breast.
5. Never use the breast pump while you are pregnant, as pumping can induce labor.
6. Inspect all appropriate pump components before each use.
7. Clean and sanitize all parts that come in contact with your breast and breastmilk prior to first use.
8. Wash all parts that come in contact with your breast and breastmilk after every use.
9. For hygienic reasons, this product is intended for use by a single user.
10. Do not allow children or pets to play with the motor unit, the adapter or accessories.
11. Remove the batteries if they are not going to be used for an extended period of time.
12. Although the breast pump is compliant with applicable EMC directives, it may still be susceptible to excessive emissions and/or may interfere with other equipment. A consequence can be that the breast pump turns off or goes into error mode. To prevent interference, keep other electric equipment away from the breast pump during expressing.
13. Do not use the device in the MR environment.


Power adapter use attention points

Use only the power adapter that comes with the product.

Make sure the voltage of the power adapter is compatible with the power source.

This product should never be left unattended when plugged into an electrical outlet.

Always unplug the breast pump immediately after use.

① 1×Seal cover

② 1×Connection tube

③ 1×Elastic diaphragm

④ 1×Seal seat

⑤ 4×Seal hoops

⑥ 1×Breastshield

⑦ 1×Pump body

⑧ 2×Valves

⑨ 4×White membranes

⑩ 2×Breastmilk bottles

⑪ 1×Bottle stand

⑫ 1×Motor unit

⑬ 2×Plugs

⑭ 1×Dome cap

⑮ 1×Nipple

⑯ 1×Screw ring

⑰ 2×Lids

⑱ 1×Power Adaptor

Installation Instructions

6

Product Assembly

Note: Make sure you have cleaned and optionally disinfected the appropriate parts of the breast pump.



After cleaning, follow these steps to assemble your collection units:

Step 1:


Insert the breastshield into the funnel-shaped section of the pump body.

Two breastshield sizes are available: 24mm and 27mm.

You can find the size on the breastshield itself.

24mm


27mm

Step 2:


Please install the white membrane on valve smoothly.

Make sure the white membrane won't curl up.


Step 3:

Insert the valve into the pump body from underneath. Push the valve in as far as possible.

Step 4:

Screw the pump body clockwise onto the bottle until it is securely fixed.

Step 5:

Put the Elastic diaphragm into the seal seat,make sure the daphragm edge and seal seat edge close fitting.Then seal by the seal cover make it as a sealed isolation assembly.If there is no "O" - shaped sealing ring on the plug at the bottom of the sealing seat ,please set the seal hoop into the groove in the plug position.


Figure 5-1


Figure 5-2


Figure 5-3


Figure 5-4

Step 6:

Make one side of connection tube without plug insert into seal cover and make another side of connection tube insert into motor unit. Then insert the whole isolation assembly into connector.


Figure 6-1


Figure 6-2


Figure 6-3



Figure 6-4

Battery Installation

This product can use 4 AA alkaline batteries as power supply. Please install the battery according to the direction of the battery polarity shown at the bottom of the unit.

Do not reverse polarity.

Do not mix old and new batteries.

Do not mix alkaline, standard(carbon-zinc) or rechargeable batteries.


Each set of batteries will provide approximately 1.5 hours of pumping time. For battery operation, the motor unit LCD will display the battery symbol .

If the battery symbol flashes, which means the battery remaining is not enough for breast pump normal work, please replace the batteries.

The Power Adapter


Although the design of the product can be used in the battery inside the unit power adapter, meanwhile, we recommend that remove the battery before use the power adapter. If you need to use power adapter when you using battery, please shut down the breast pump then access external power supply and turn on the product.

When the external power supply is used, the battery symbol will not show on the screen. At the same time the external power symbol  will light up.


Please don't insert or pull out the power adapter with wet hand.
Please don't plug the power adapter with the power supply for a long time.

Display information and button function


- ① The external power supply indicator
- ② Pumping level indicator
- ③ Continuity work time, "minute" in front, "second" behind
- ④ Stimulation mode
- ⑤ Expression mode
- ⑥ Battery indicator
- ⑦ Level increasing direction

Relax in a comfortable chair (you may wish to use cushions to support your back).

- a. Press the assembled pump body against your breast. Make sure that your nipple is centered. Hold the breastshield onto your breast with your thumb and index finger. Support your breast with the palm of your hand.

- b. Press the  button, the breast pump automatically starts up in stimulation mode as per preset level and the stimulation indicator lights up. If the breast pump used before without cut off power supply, the breast pump will express at last set suction level.
- c. Once your milk starts flowing, please press the “Stimulation/Expression” button to enter into expression mode. You can switch the stimulation mode and expression mode by press the “Stimulation/Expression” button as per your requirement. If no action is taken during stimulation mode, after two minutes, the pump will automatically change to the expression mode.
- d. Depending on your own personal comfort, you can press the “+” or “-” button to adjust the suction level. Long press the button can help you adjust the level quicker.

- e. Close the bottle with a lid after pumping. Turn off the breast pump.
Disassemble and clean the parts that contact breast and milk.



Please keep the airway clear and unobstructed during sucking, prevent air path blocked.

Use bottle holder to prevent the bottle from tipping over. Do not fill the bottle too full to prevent overfilling and spillage.

The unit can turn off the power itself about 30 minutes no operation.

Make sure that your nipple fits correctly in the breast shield. Not too tight, it is important that the nipple can move freely, while expressing breast milk. The breastshield should not be too large, because this can cause pain or less effective expression. Two breast shield sizes are available: 24mm and 27mm. There's one size breast shield come with the product.

24mm

27mm

If there is any pain or discomfort during pumping, please consider to choose a larger or smaller size.

Caution:

1. Only use drinking-quality tap or bottled water for cleaning.
2. Take apart and wash all parts that come in contact with the breast and breastmilk immediately before and after use to avoid dry up of milk residues and to prevent growth of bacteria.
3. Please disassemble all parts which require clean to make sure cleaning thoroughly.
4. Please place parts on a clean environment to avoid being polluted.
5. Never put breastpump motor unit and AC adapter in water or a sterilizer, as you can cause permanent damage to the breastpump.

Cleaning before use:

1. Take apart and wash all parts that come in contact with the breast and breastmilk.
2. Put all separated parts into pot. Fill in the pot with enough drinking-quality tap or bottled water to cover all parts.
3. Bring the water to boil. Place parts in boiling water for 5 minutes.
4. Allow water to cool and gently remove the parts from the water. Be careful and don't scald your skin.
5. Place the parts neatly on a clean paper towel or in a clean drying rack and allow them to air dry. Avoid using cloth towels to dry the parts because they can carry germs and bacteria that are harmful to your baby.

Cleaning after use:

1. Take apart and wash all parts that come in contact with the breast and breastmilk.
2. Rinse in cool water all separated parts that came in contact with breast and breastmilk in order to remove breastmilk residue.
3. Place the parts neatly on a clean paper towel or in a clean drying rack and allow them to air dry.

Abnormal phenomenon	Cause analysis	Solutions
Low or No Suction	Connection points is loose	Inspect all connection points to ensure attachment is secure
	The valve is chipped or cracked or there's holes or tears on membrane is	Replace valve and/or membrane prior to pumping
	Breastshield size is not fit	Replace a appropriate size breastshield
Shows "Err"	Power problems	Please replace batteries or power adapter.
Battery symbol flash	Voltage is too low	Please replace batteries


Correct Disposal of This Product


(Waste Electrical & Electronic Equipment)

This marking shown on the product indicates that it should not be disposed with other household waste at the end of its life. To prevent potential harm to the environment or to human health, please separate this product from other types of wastes and recycle it responsibly. When disposing this type of product, contact the retailer where product was purchased or contact your local government office for details regarding how this item can be disposed in an environmentally safe recycling center. Business users should contact their supplier and check the terms and conditions of the purchasing agreement. This product should not be mixed with other commercial wastes for disposal. This product is free of hazardous materials.

1. Avoid dropping, slamming, or throwing the unit.

2. Avoid extreme temperatures. Do not expose unit directly under sunshine.


3. When cleaning the unit, use a soft fabric and lightly wipe with mild detergent.


Use a damp cloth to remove dirt and excess detergent.


4. Do not use petrol, thinners or similar solvents.

5. Remove batteries when not in operation for an extended period of time.

6. Do not disassemble product.

Product Description	Electric Breast Pump	
Model	Smilky	
Intermittent Vacuum	Stimulate mode	Approx. -60~-130mmHg
	Expression mode	Approx. -60~-250mmHg
Vacuum Levels	Stimulate mode	10 Levels
	Expression mode	10 Levels
Cycles Per Minute	Stimulate mode	100 C.P.M.
	Expression mode	29~73 C.P.M.
Breastshield Sizes	24mm	
Power Source	4 AA batteries or Medical AC Adapter(DC6.0V, 1000mA)	
Battery Life	Approximately 1.5 hours,different types of batteries may affect battery life	
Additional Function	Automatic Power-Off	
Unit Weight	Approx. 265g(9.35oz)(excluding battery)	
Unit Dimensions	129 x 129 x 55mm (L x W x H)	
Operating Environment	Temperature	5°C ~ 40°C(41°F~104°F)
	Humidity	15~93%RH
	Pressure	700hPa~1060hPa
Storage Environment	Temperature	-20°C~-55°C(-4°F~131°F)
	Humidity	≤93%RH
Expected Service Life	Your appliance has been designed and developed with the greatest possible care and has an expected service life of 400 hours.	

Ingress Protection Rating	IP21
Classification	Internal Powered Equipment, Type BF 

Specifications are subject to change without notice.

This appliance conforms to the following standards:

EN 60601-1-11 Medical electrical equipment —Part 1-11: General requirements for basic safety and essential performance – Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment and complies with the requirements of EN 60601-1-2(EMC), IEC/EN60601-1(Safety) standards. And the manufacturer is ISO 13485 certified.

We warrant this product to be free from manufacturing defects for a period of one year for pump mechanism (90 days for parts and accessories) from the original date of consumer purchase. We gladly offer a free exchange and lifetime repair due to manufacturing defects. This warranty does not apply to batteries and any product which has been subjected to misuse, abuse or alteration. Please contact local retailer for details.

Contact Information

The lay operator or lay responsible or ganization should contact the manufacturer or the representative of manufacturer.

- for assistance,if needed,in setting up,using or maintaining the product,or
- to report unexpected operation or events.

The device satisfies the EMC requirements of the international standard IEC 60601-1-2. The requirements are satisfied under the conditions described in the table below. The device is an electrical medical product and is subject to special precautionary measures with regard to EMC which must be published in the instructions for use. Portable and mobile HF communications equipment can affect the device. Use of the unit in conjunction with non-approved accessories can affect the device negatively and alter the electromagnetic compatibility. The device should not be used directly adjacent to or between other electrical equipment.

Table 1

Guidance and declaration of manufacturer-electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment -guidance
Radiated emission CISPR 11	Group 1, class B.	The device uses RF energy only for its internal function. Therefore, its emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Conducted emission CISPR 11	Group 1, class B.	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Table 2

Guidance and declaration of manufacturer-electromagnetic immunity (For home healthcare environment)			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment -guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrostatic transient/burst IEC 61000-4-4	± 2 kV, 100kHz, for AC power port	± 2 kV, 100kHz, for AC power port	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±0.5kV, ±1kV (differential mode)	±0.5kV, ±1kV (differential mode)	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225° , 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycle	0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225° , 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m; 50Hz or 60Hz	30 A/m; 50Hz or 60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Table 3

Guidance and declaration of manufacturer-electromagnetic immunity (For home healthcare environment)			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment -guidance
Conducted RF IEC 61000-4-6	3V for 0.15-80MHz; 6V in ISM and amateur radio bands between 0.15-80MHz	3V for 0.15-80MHz; 6V in ISM and amateur radio bands between 0.15-80MHz	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
	Radiated RF IEC 61000-4-3	385MHz, 27V/m	
	450MHz, 28V/m	450MHz, 28V/m	Recommended separation distance $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ 80 MHz to 800 MHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$ 800 MHz to 2.7 GHz
	710MHz, 745MHz, 780MHz 9V/m	710MHz, 745MHz, 780MHz 9V/m	where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).
	810MHz, 870MHz, 930MHz 28V/m	810MHz, 870MHz, 930MHz 28V/m	
	1720MHz, 1845MHz, 1970MHz 28V/m	1720MHz, 1845MHz, 1970MHz 28V/m	Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.
	2450MHz, 28V/m	2450MHz, 28V/m	Interference may occur in the vicinity of equipment marked with the following symbol:
	5240MHz, 5500MHz, 5785MHz 9V/m	5240MHz, 5500MHz, 5785MHz 9V/m	

Table 4

Recommended separation distances between portable and mobile RF communications equipment and the device (For home healthcare environment)		
The device is intended for use in an electromagnetic environment in which radiated therefore disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.		
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m	
	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1}\right]\sqrt{P}$	800 MHz to 2.7 GHz $d = \left[\frac{7}{E_1}\right]\sqrt{P}$
0.01	0.12	0.23
0.1	0.38	0.73
1	1.2	2.3
10	3.8	7.3
100	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance <i>d</i> in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where <i>P</i> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.		
NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.		
NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		

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