Nebulizer

health enjoyment

LUCKY

Nebulizer - COMPRESSOR TYPE

LUCKY has adorable, cute, child-friendly design. Operation of nebulizer is very easy with just one buton. Suitable for adult and pediatric patients.

- One button easy operation
 - 2 Effective medication delivery
 - O Adorable and friendly design
 - I First choice, ideal for all children





Specifications

Intended Use	Aerosolization of drugs for inhalation by the patient
Environment of Use	Home, Hospital, Sub -acute Institutions
Patient Population	Pediatric
Medication Capacity	6 ml
MMAD:	3µm
Average Nebulization Rate:	Min 0.32ml / Min
Noise Level:	B elow 58 dB(A)
Compressor Pressure Range:	35 - 50Psi
Max. Operating Pressure :	9-16 Psi (62kpa to 110 kpa)
Flow Rate:	5~8 lpm
Mode of operation	Continuous up to 4 hours
Contraindications	None
Protection against electric shock	- classified as Type II
Type BF applied parts:	Mouthpiece, masks
IP Class	IP21
Materials in patient contact	Polypropylene
Nebulizer components cleanable	Yes
Operating conditions	5°C - 40°C, 10% - 95% RH
Storage conditions	(-25°C) – (+60°C), 10% - 95% RH
Power:	230V AC, 50Hz
Power Consumption:	Below 65W
Fuse:	2A
Dimensions (mm) (WxDxH)	190 x 220 x 150
Weight (kg)	1.8 kg
Standard Accessories:	Compressor Nebulizer Machine : 1pcs Adult Mask: 1pcs Child Mask: 1pcs Air Filter:1pcs Durable Air Tube:1pcs Medical Cup 6ml:1pcs
	Mouth Piece: 1pcs Nosepiece: 1pcs

Descriptions



- 1. Nebulizer machine
- 2. Nebulizer cup
- 3. Air tube
- 4. Mouthpiece
- 5. Air filters
- 6. Adult mask
- 7. Pediatric mask
- 8. Nosepiece

www.trimpeks.com