

## STEAMY

## **Nebulizer - COMPRESSOR TYPE**

"Steamy" is a simple and easy to use nebulizer with a powerful compressor. It's designed to deliver the prescribed medication solution to treat patient respiratory disorders, such as asthma, allergies and bronchitis. The nebulizer converts the medication solution into an aerosolized mist which is inhaled by the patient through the mouthpiece or masks.

- Powerful compressor
- 2 6cc. nebulizer cup
- 3 Particle size 3 µm
- 4 Adult and child masks
- 6 Nebulizer kit



## **Features**

T	
Intended Use	Aerosolization of drugs for inhalation by the patient
<b>Environment of Use</b>	Home, Hospital, Sub-acute Institutions
Patient Population	Pediatric and adult
Medication Capacity	6 ml
MMAD:	3μm
Average Nebulization Rate:	Min 0.32ml / Min
Noise Level:	Below 58 dB(A)
<b>Compressor Pressure Range:</b>	35 - 50Psi
Max. Operating Pressure:	9-16 Psi ( 62kpa to 110 kpa)
Flow Rate:	5∼8 lpm
Mode of operation	Continuous up to 4 hours
Contraindications	None
Protection against electric shock	- classified as Type II
Type BF applied parts:	Mouthpiece, masks
IP Class	IP21
Materials in patient contact	Polypropylene
Nebulizer components cleanable	Yes
Operating conditions	5°C - 40°C, 10% - 95% RH
Storage conditions	(-25°C) – (+60°C), 10% - 95% RH
Power:	230V AC, 50Hz
<b>Power Consumption:</b>	Below 65W
Fuse:	2A
Dimensions (mm) (WxDxH)	280 x 165 x 90
Weight (kg)	1.80 kg
Standard Accessories:	Compressor Nebulizer Machine: 1pcs Adult Mask: 1pcs Child Mask: 1pcs Air Filter: 1pcs Durable Air Tube: 1pcs
	Medical Cup 6ml:1pcs Mouth Piece: 1pcs Nosepiece: 1pcs





- 1. Nebulizer machine
- 2. Nebulizer cup
- 3. Air tube
- 4. Mouthpiece
- 5. Air filters
- 6. Adult mask
- 7. Pediatric mask
- 8. Nosepiece